



	Instant (< 10 minutes)	Short term (0-1 day)	Medium term (1-3 days)	Long term (0-7 days)
Clinical Applications	<ul style="list-style-type: none"> Urine sample collection Check of passage in urethra Bladder tension relief Bladder lavage 	<ul style="list-style-type: none"> Urine sample collection Check of passage in urethra Bladder tension relief Bladder lavage 	<ul style="list-style-type: none"> Urine sample collection Check of passage in urethra Bladder tension relief Bladder lavage Catheterization for drainage purposes during surgery or hospitalization 	<ul style="list-style-type: none"> Urine sample collection Bladder tension relief Bladder lavage Catheterization for drainage purposes during surgery or hospitalization
Recommended Catheters	<p>BUSTER HC Urinary Catheter</p> <ul style="list-style-type: none"> ✓ No need for gel ✓ Soft and comfortable ✓ Colour-coded hub and label for easy size recognition ✓ Luer lock on hub allows easy connection to syringe 	<p>BUSTER Male Dog Catheter</p> <ul style="list-style-type: none"> ✓ More rigid catheter ✓ Luer lock on hub allows easy connection to syringe or urine bag 	<p>BUSTER NC Urinary Catheter</p> <ul style="list-style-type: none"> ✓ Soft and comfortable ✓ Colour-coded hub and label for easy size recognition ✓ Luer lock on hub allows easy connection to syringe or urine bag 	<p>BUSTER Foley Catheter</p> <ul style="list-style-type: none"> ✓ Very soft silicone catheter with high biocompatibility ✓ Colour-coded hub and label for easy size recognition ✓ Urine drainage port fits KRUUSE urine bag ✓ Radiopaque for x-ray visibility